

Material description: PETROLEUM SPIRIT 80/100 pure
CAS No: 8032-32-4

<i>Supplied as:</i>	Product code	Description
	CHE2814	PETROLEUM SPIRIT 80/100 pure 1L.
	CHE2816	PETROLEUM SPIRIT 80/100 pure 2.5L.
	CHE2818	PETROLEUM SPIRIT 80/100 pure 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Boiling range; 90%	°C	80	100		
Density @ 15°C	g/ml	0.690	0.731		

Revision date: 04/09/2024

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1304/MS/20240904