



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: PHENOL A.R.

CAS No: 108-95-2

Supplied as: Product code
CHE2822

Description
PHENOL A.R. 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Melting point	°C	40.5			
Boiling point	°C			181	
Moisture	%		0.5		
Non-volatiles	%		0.05		

Revision date: 24/09/1991

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1306/MS/19910924