



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: PHENOL RED SOLUTION

CAS No: Mixture

Supplied as: Product code Description
CHE2830 PHENOL RED SOLUTION 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Yellow to Red					pH 6.8-8.4

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1309/MS/19890101