



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: PHLOROGLUCINOL pure
CAS No: 108-73-6

Supplied as: Product code Description
CHE2860 PHLOROGLUCINOL pure 10g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	95			
Melting range; dried	°C	218	221		

Revision date: 08/12/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1944/MS/20201208