



Material description: PHTHALIC ANHYDRIDE pure
CAS No: 85-44-9

Supplied as:	Product code	Description
	CHE2872	PHTHALIC ANHYDRIDE pure 500g.
	CHE4044	PHTHALIC ANHYDRIDE pure 1Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.0			
Melting range	°C	132	134		

Revision date: 01/01/1989 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1329/MS/19890101