



Material description: POTASSIUM CHLORIDE 4M (pH electrode filling soln)

CAS No: Mixture

Supplied as: Product code Description
CHE2910 POTASSIUM CHLORIDE 4M (pH electrode filling soln) 100ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Density @20°C	g/ml	1.172	1.176	1.174	

Revision date: 09/07/1998

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2859/MS/19980709