

Material description: POTASSIUM HYDROXIDE PELLETS A.R.

CAS No: 1310-58-3

<i>Supplied as:</i>	Product code	Description
	CHE3000	POTASSIUM HYDROXIDE PELLETS A.R. 500g.
	CHE3002	POTASSIUM HYDROXIDE PELLETS A.R. 1Kg.
	CHE3976	POTASSIUM HYDROXIDE PELLETS A.R. 100g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	85			
Carbonate	%		2.0		
Chloride	%		0.020		
Phosphate	%		0.0100		
Sulphate	%		0.020		
Iron	%		0.001		
Aluminium	%		0.00002		
Heavy Metals; as Pb	%		0.0010		
Identity					Complies with Ph.Eur/BP
Appearance of aqueous solution					Complies with Ph.Eur/BP

Revision date: 28/11/2019

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0539/MS/20191128