

Material description: POTASSIUM HYDROXIDE PELLETS pure
CAS No: 1310-58-3

<i>Supplied as:</i>	Product code	Description
	CHE3004	POTASSIUM HYDROXIDE PELLETS pure 500g.
	CHE3006	POTASSIUM HYDROXIDE PELLETS pure 1Kg.
	CHE3008	POTASSIUM HYDROXIDE PELLETS pure 2Kg.
	CHE3009	POTASSIUM HYDROXIDE PELLETS pure 25Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	85			
Carbonate	%		2.0		
Chloride	%		0.02		
Sulphate	%		0.02		

Revision date: 28/11/2019 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0538/MS/20191128