



Material description: POTASSIUM HYDROXIDE SOLUTION 70% w/v (50% w/w)
CAS No: 1310-58-3

Supplied as: Product code Description
CHE3012 POTASSIUM HYDROXIDE SOLUTION 70% w/v (50% w/w) 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	% w/w	49	50		
Wt. per ml.	g			1.50	
Iron	ppm		5.0		

Revision date: 27/01/2021 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1349/MS/20210127