

Material description: POTASSIUM IODATE A.C.S.

CAS No: 7758-05-6

Supplied as: Product code Description
CHE3018 POTASSIUM IODATE A.C.S. 500g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5	100.4		
pH; 5% solution	pH	5.0	9.0		
Iodide	%		0.001		
Other halides	%		0.02		
Sulphates	%		0.01		
Iron	%		0.001		
Sodium	%		0.01		
Moisture	%		0.05		
Insoluble matter	%		0.005		

Revision date: 10/04/2006

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2005/MS/20060410