

Material description: POTASSIUM IODATE (V) pure  
CAS No: 7758-05-6

<i>Supplied as:</i>	Product code	Description
	CHE3020	POTASSIUM IODATE (V) pure 100g.
	CHE3022	POTASSIUM IODATE (V) pure 250g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Chlorate	%		0.05		
Iodide	%		0.003		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0540/MS/19890101