

Material description: POTASSIUM IODIDE A.R.

CAS No: 7681-11-0

<i>Supplied as:</i>	Product code	Description
	CHE3028	POTASSIUM IODIDE A.R. 500g.
	CHE3030	POTASSIUM IODIDE A.R. 100g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity; dried	%	99.00			
Sulphate					Complies BP
Iron	ppm		20		
Heavy Metals; as Pb	ppm		10		
Loss on drying	%		1.0		
Iodates					Complies BP
Thiosulphates, Barium					Complies BP
Nitrate, Nitrite, Ammonia					Complies BP
Solubility					Complies BP
Identity					Complies BP
Appearance of aqueous solution					Complies BP

Revision date: 21/08/2014

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0542/MS/20140821