



Material description: POTASSIUM IODIDE SOLUTION 10% w/v

CAS No: Mixture

Supplied as: Product code

CHE3042

Description

POTASSIUM IODIDE SOLUTION 10% w/v 500ml.

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|---------------|-------|---------------|---------------|---------------|-------|
| Density @20°C | g/ml | 1.064 | 1.074 | 1.069 | |

Revision date: 18/09/2012

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1859/MS/20120918