

Material description: POTASSIUM OXALATE H2O pure  
CAS No: 6487-48-5

<i>Supplied as:</i>	Product code	Description
	CHE3058	POTASSIUM OXALATE H2O pure 250g.
	CHE3062	POTASSIUM OXALATE H2O pure 1Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Chloride	%		0.01		
Iron	%		0.005		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1352/MS/19890101