

Material description: POTASSIUM PERMANGANATE A.R.

CAS No: 7722-64-7

|                     |              |                                   |
|---------------------|--------------|-----------------------------------|
| <i>Supplied as:</i> | Product code | Description                       |
|                     | CHE3072      | POTASSIUM PERMANGANATE A.R. 500g. |
|                     | CHE3074      | POTASSIUM PERMANGANATE A.R. 250g. |

## Specification

| Parameter        | Units | Minimum value | Maximum value | Typical value | Notes |
|------------------|-------|---------------|---------------|---------------|-------|
| Purity           | %     | 99.5          |               |               |       |
| Chloride         | %     |               | 0.005         |               |       |
| Nitrate          | %     |               | 0.04          |               |       |
| Sulphate         | %     |               | 0.01          |               |       |
| Arsenic          | %     |               | 0.0001        |               |       |
| Insoluble matter | %     |               | 0.05          |               |       |

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0549/MS/19890101