

Material description: POTASSIUM PERMANGANATE A.R.

CAS No: 7722-64-7

<i>Supplied as:</i>	Product code	Description
	CHE3072	POTASSIUM PERMANGANATE A.R. 500g.
	CHE3074	POTASSIUM PERMANGANATE A.R. 250g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Chloride	%		0.005		
Nitrate	%		0.04		
Sulphate	%		0.01		
Arsenic	%		0.0001		
Insoluble matter	%		0.05		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0549/MS/19890101