



Material description: POTASSIUM PERMANGANATE 0.02M (N/10)

CAS No: Mixture

Supplied as: Product code

CHE3078

Description

POTASSIUM PERMANGANATE 0.02M (N/10) 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	0.01998	0.02002		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1783/MS/19890101