

Material description: PROPAN-1-OL pure  
CAS No: 71-23-8

<i>Supplied as:</i>	Product code	Description
	CHE3104	PROPAN-1-OL pure 2.5L.
	CHE3106	PROPAN-1-OL pure 250ml.
	CHE3108	PROPAN-1-OL pure 500ml.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Wt. per ml.	g	0.802	0.807		
Acidity; acetic	%		0.003		
Water	% w/w		0.3		
Colour	Hazen		10		

Revision date: 04/01/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1362/MS/20210104