

Material description: PROPAN-2-OL HPLC
CAS No: 67-63-0

Supplied as: Product code Description
CHE3114 PROPAN-2-OL HPLC 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.5			
Water	% w/w		0.10		
Residue on evaporation	ppm		2.00		
Absorbance @ 220nm	AU		0.200		
Absorbance @ 230nm	AU		0.100		
Absorbance @ 240nm	AU		0.050		
Absorbance @ 250nm	AU		0.010		
Absorbance @ 260nm	AU		0.005		
Absorbance @ 270nm	AU		0.005		
Wavelength @ 1AU	nm		205		

Revision date: 23/01/2020 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2232/MS/20200123