



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: PROPAN-2-OL pure
CAS No: 67-63-0

| | | |
|---------------------|--------------|-------------------------|
| <i>Supplied as:</i> | Product code | Description |
| | CHE3116 | PROPAN-2-OL pure 2.5L. |
| | CHE3118 | PROPAN-2-OL pure 500ml. |
| | CHE3120 | PROPAN-2-OL pure 25L. |

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|---------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 99.5 | | | |
| Boiling range | °C | 81 | 83 | | |
| Wt. per ml. | g | 0.784 | 0.787 | | |

Revision date: 04/02/1994

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0557/MS/19940204