

Material description: SAND PURIFIED (dried white silica sand)

CAS No: 7631-86-9

<i>Supplied as:</i>	Product code	Description
	CHE3184	SAND PURIFIED (dried white silica sand) 1Kg.
	CHE3186	SAND PURIFIED (dried white silica sand) 5Kg.
	CHE3188	SAND PURIFIED (dried white silica sand) 25Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Silica	%			99.7	
Aluminium Oxide	%			0.15	
Ferric Oxide	%			0.01	
Sodium Oxide	%		0.03		
Potassium Oxide	%			0.050	
Calcium Oxide	%			0.050	
Titanium Dioxide	%			0.020	
Acid washed					

Revision date: 14/02/2024

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1788/MS/20240214