



Material description: SILVER NITRATE 0.02M (N/50)

CAS No: Mixture

<i>Supplied as:</i>	Product code	Description
	CHE3238	SILVER NITRATE 0.02M (N/50) 2.5L.
	CHE3240	SILVER NITRATE 0.02M (N/50) 1L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	0.01998	0.02002		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2035/MS/19890101