

Material description: SILVER NITRATE SOLUTION 2% w/v

CAS No: Mixture

Supplied as: Product code Description
CHE3250 SILVER NITRATE SOLUTION 2% w/v 500ml.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Silver nitrate	g/L			20	
Nitric acid	mol/L			0.003	
Density @20°C	g/ml	1.022	1.030		

Revision date: 03/09/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3088/MS/20200903