



Material description: SODA LIME SELF-INDICATING (3-5mm)

CAS No: 8006-28-8

Supplied as: Product code Description
CHE3256 SODA LIME SELF-INDICATING (3-5mm) 500g.

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|--------------------------|-------|---------------|---------------|---------------|-------------|
| Colour change; anhydrous | | | | | Off white |
| Colour change; hydrated | | | | | Violet blue |

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2050/MS/19890101