

Material description: SODIUM CHLORIDE A.R.

CAS No: 7647-14-5

<i>Supplied as:</i>	Product code	Description
	CHE3314	SODIUM CHLORIDE A.R. 3Kg.
	CHE3316	SODIUM CHLORIDE A.R. 500g.
	CHE3318	SODIUM CHLORIDE A.R. 25Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99	101		
Moisture	%		0.5		
Bromide	%		0.01		
Phosphate	%		0.0025		
Aluminium	%		0.00002		
Arsenic	%		0.0001		
Iron	%		0.0002		
Lead	%			0.0005	
Nitrites					Complies BP

Revision date: 07/05/2021

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0589/MS/20210507