

Material description: SODIUM DIHYDROGEN ORTHOPHOSPHATE 2H<sub>2</sub>O pure  
CAS No: 13472-35-0

<i>Supplied as:</i>	Product code	Description
	CHE3330	SODIUM DIHYDROGEN ORTHOPHOSPHATE 2H <sub>2</sub> O pure 250g.
	CHE3332	SODIUM DIHYDROGEN ORTHOPHOSPHATE 2H <sub>2</sub> O pure 500g.
	CHE3334	SODIUM DIHYDROGEN ORTHOPHOSPHATE 2H <sub>2</sub> O pure 3Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Arsenic	%		0.0002		
Sulphate	%		0.15		
Iron	%		0.002		

Revision date: 01/01/1989      Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1947/MS/19890101