

Material description: SODIUM HYDROGEN CARBONATE A.R.

CAS No: 144-55-8

<i>Supplied as:</i>	Product code	Description
	CHE3350	SODIUM HYDROGEN CARBONATE A.R. 500g.
	CHE3352	SODIUM HYDROGEN CARBONATE A.R. 3Kg.
	CHE3354	SODIUM HYDROGEN CARBONATE A.R. 1Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.8			
Chloride	%		0.001		
Sulphate	%		0.003		
Arsenic	%		0.00004		
Iron	%		0.0005		
Lead	%		0.0002		
Potassium	%		0.01		
Insoluble matter	%		0.02		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0601/MS/19890101