

Material description: SODIUM HYDROGEN CARBONATE A.R.

CAS No: 144-55-8

<i>Supplied as:</i>	Product code	Description
	CHE3350	SODIUM HYDROGEN CARBONATE A.R. 500g.
	CHE3352	SODIUM HYDROGEN CARBONATE A.R. 3Kg.
	CHE3354	SODIUM HYDROGEN CARBONATE A.R. 1Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.8			
Chloride	%		0.001		
Sulphate	%		0.003		
Arsenic	%		0.00004		
Iron	%		0.0005		
Lead	%		0.0002		
Potassium	%		0.01		
Insoluble matter	%		0.02		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0601/MS/19890101