

Material description: SODIUM HYDROGEN CARBONATE A.R.

CAS No: 144-55-8

| <i>Supplied as:</i> | Product code | Description |
|---------------------|--------------|--------------------------------------|
| | CHE3350 | SODIUM HYDROGEN CARBONATE A.R. 500g. |
| | CHE3352 | SODIUM HYDROGEN CARBONATE A.R. 3Kg. |
| | CHE3354 | SODIUM HYDROGEN CARBONATE A.R. 1Kg. |

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|------------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 99.8 | | | |
| Chloride | % | | 0.001 | | |
| Sulphate | % | | 0.003 | | |
| Arsenic | % | | 0.00004 | | |
| Iron | % | | 0.0005 | | |
| Lead | % | | 0.0002 | | |
| Potassium | % | | 0.01 | | |
| Insoluble matter | % | | 0.02 | | |

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0601/MS/19890101