

Material description: SODIUM HYDROGEN SULPHATE ANHYDROUS pure
CAS No: 7681-38-1

<i>Supplied as:</i>	Product code	Description
	CHE3362	SODIUM HYDROGEN SULPHATE ANHYDROUS pure 500g.
	CHE3364	SODIUM HYDROGEN SULPHATE ANHYDROUS pure 1Kg.
	CHE3366	SODIUM HYDROGEN SULPHATE ANHYDROUS pure 3Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	92			
Chloride	%			0.05	
Iron	%			0.02	
Moisture	%			1	
Insoluble matter	%			0.02	

Revision date: 18/10/1993 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0605/MS/19931018