



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: SODIUM HYDROXIDE 0.111M (N/9)

CAS No: 1310-73-2

| <i>Supplied as:</i> | Product code | Description |
|---------------------|--------------|-------------------------------------|
| | CHE3406 | SODIUM HYDROXIDE 0.111M (N/9) 5L. |
| | CHE3408 | SODIUM HYDROXIDE 0.111M (N/9) 2.5L. |

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|-----------|-------|---------------|---------------|---------------|-------|
| Strength | mol/L | 0.11089 | 0.11111 | | |

Revision date: 01/01/1989

Issue date: 12/05/2022

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1796/MS/19890101