

Material description: SODIUM HYDROXIDE PEARL pure  
CAS No: 1310-73-2

<i>Supplied as:</i>	Product code	Description
	CHE3410	SODIUM HYDROXIDE PEARL pure 500g.
	CHE3412	SODIUM HYDROXIDE PEARL pure 25Kg.
	CHE3414	SODIUM HYDROXIDE PEARL pure 1Kg.
	CHE4102	SODIUM HYDROXIDE PEARL pure (SUPPLIERS PACKAGING) 25Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98			
Chloride	%		0.02		
Sulphate	%		0.1		
Iron	%		0.002		

Revision date: 03/01/1992      Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1419/MS/19920103