

Material description: SODIUM HYDROXIDE PELLETS A.R.

CAS No: 1310-73-2

<i>Supplied as:</i>	Product code	Description
	CHE3416	SODIUM HYDROXIDE PELLETS A.R. 500g.
	CHE3418	SODIUM HYDROXIDE PELLETS A.R. 3Kg.
	CHE3420	SODIUM HYDROXIDE PELLETS A.R. 1Kg.
	CHE3972	SODIUM HYDROXIDE PELLETS A.R. 100g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	97.5			
Carbonate	%		1.0		
Chloride	%		0.01		
Sulphate	%		0.01		
Lead	ppm		20.0000		
Iron	%		0.001		

Revision date: 04/01/2011

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0608/MS/20110104