

Material description: SODIUM HYDROXIDE PELLETS A.R.

CAS No: 1310-73-2

<i>Supplied as:</i>	Product code	Description
	CHE3416	SODIUM HYDROXIDE PELLETS A.R. 500g.
	CHE3418	SODIUM HYDROXIDE PELLETS A.R. 3Kg.
	CHE3420	SODIUM HYDROXIDE PELLETS A.R. 1Kg.
	CHE3972	SODIUM HYDROXIDE PELLETS A.R. 100g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	97.5			
Carbonate	%		1.0		
Chloride	%		0.01		
Sulphate	%		0.01		
Lead	ppm		20.0000		
Iron	%		0.001		

Revision date: 04/01/2011

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0608/MS/20110104