



Material description: SODIUM HYDROXIDE SOLUTION 70% w/v (48% w/w)  
CAS No: Mixture

Supplied as: Product code Description  
CHE3450 SODIUM HYDROXIDE SOLUTION 70% w/v (48% w/w) 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	% w/v	68.0	71.5		
Wt. per ml.	g			1.49	

Revision date: 10/07/1995 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1418/MS/19950710