

Material description: PURECHLOR 15 (SODIUM HYPOCHLORITE SOLUTION)

CAS No: 7681-52-9

<i>Supplied as:</i>	Product code	Description
	CHE3454	PURECHLOR 15 (SODIUM HYPOCHLORITE SOLUTION) 2.5L.
	CHE3456	PURECHLOR 15 (SODIUM HYPOCHLORITE SOLUTION) 5L.
	CHE3458	PURECHLOR 15 (SODIUM HYPOCHLORITE SOLUTION) 25L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Wt. per ml.	g			1.27	
Available chlorine	% w/w	14			
Sodium hydroxide	% w/w		1		
Sodium carbonate	% w/w		1.5		
Iron	ppm		2		

Revision date: 19/05/2009

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0609/MS/20090519