

Material description: SODIUM METABISULPHITE pure
CAS No: 7681-57-4

<i>Supplied as:</i>	Product code	Description
	CHE3474	SODIUM METABISULPHITE pure 500g.
	CHE3476	SODIUM METABISULPHITE pure 3Kg.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	95.5			
Iron	%		0.002		
Lead	%		0.001		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0610/MS/19890101