

Material description: SODIUM METAL (in paraffin liquid)

CAS No: 7440-23-5

<i>Supplied as:</i>	Product code	Description
	CHE3478	SODIUM METAL (in paraffin liquid) 25g.
	CHE3480	SODIUM METAL (in paraffin liquid) 500g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99.8			
Chloride	%		0.01		
Calcium	%		0.075		
Iron	%		0.001		
Potassium	%		0.03		

Revision date: 30/06/1994

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0612/MS/19940630