

Material description: SODIUM NITRITE pure
CAS No: 7632-00-0

| <i>Supplied as:</i> | Product code | Description |
|---------------------|--------------|---------------------------|
| | CHE3498 | SODIUM NITRITE pure 250g. |
| | CHE3500 | SODIUM NITRITE pure 500g. |
| | CHE3502 | SODIUM NITRITE pure 3Kg. |

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|------------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 98 | | | |
| Insoluble matter | % | | 0.01 | | |

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0615/MS/19890101