

Material description: SODIUM SULPHATE ANHYDROUS A.R.

CAS No: 7757-82-6

| <i>Supplied as:</i> | Product code | Description |
|---------------------|--------------|--------------------------------------|
| | CHE3538 | SODIUM SULPHATE ANHYDROUS A.R. 500g. |
| | CHE3540 | SODIUM SULPHATE ANHYDROUS A.R. 1Kg. |
| | CHE3542 | SODIUM SULPHATE ANHYDROUS A.R. 3Kg. |
| | CHE3970 | SODIUM SULPHATE ANHYDROUS A.R. 100g. |

Specification

| Parameter | Units | Minimum value | Maximum value | Typical value | Notes |
|------------------|-------|---------------|---------------|---------------|-------|
| Purity | % | 99.5 | | | |
| Chloride | % | | 0.18 | | |
| Insoluble matter | % | | 0.01 | | |
| pH; 1% solution | pH | 6.0 | 8.0 | | |

Revision date: 24/09/2020

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0623/MS/20200924