

Material description: SODIUM SULPHATE ANHYDROUS pure

CAS No: 7757-82-6

<i>Supplied as:</i>	Product code	Description
	CHE0622-20	SODIUM SULPHATE ANHYDROUS pure 10Kg.
	CHE3544	SODIUM SULPHATE ANHYDROUS pure 500g.
	CHE3546	SODIUM SULPHATE ANHYDROUS pure 3Kg.
	CHE3548	SODIUM SULPHATE ANHYDROUS pure 1Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Chloride	%		0.2		

Revision date: 22/10/2015

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0622/MS/20151022