

Material description: SODIUM TAUROGLYCOCHOLATE (BILE SALTS)
CAS No: 11006-55-6

Supplied as: Product code Description
CHE3556 SODIUM TAUROGLYCOCHOLATE (BILE SALTS) 25g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Identity					Complies
Bile Acid Content	%	45.0			
Loss on drying	%		5.00		
pH ;2% solution	pH	6.00	8.00		

Revision date: 01/09/2022

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 3277/MS/20220901