

Material description: SODIUM THIOSULPHATE 5H<sub>2</sub>O A.R.

CAS No: 10102-17-7

<i>Supplied as:</i>	Product code	Description
	CHE3566	SODIUM THIOSULPHATE 5H <sub>2</sub> O A.R. 500g.
	CHE3568	SODIUM THIOSULPHATE 5H <sub>2</sub> O A.R. 3Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	98.5	101.0		
pH ;10% solution	pH	6.5	9.0		
Alkali metals; as Ca	%		0.2		

Revision date: 11/10/2023

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0632/MS/20231011