

Material description: TIN (II) CHLORIDE 2H<sub>2</sub>O pure  
CAS No: 10025-69-1

<i>Supplied as:</i>	Product code	Description
	CHE3608	TIN (II) CHLORIDE 2H <sub>2</sub> O pure 100g.
	CHE3610	TIN (II) CHLORIDE 2H <sub>2</sub> O pure 250g.
	CHE3612	TIN (II) CHLORIDE 2H <sub>2</sub> O pure 500g.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	97.5			
Iron	%		0.005		
Lead	%		0.003		

Revision date: 29/08/1991

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0635/MS/19910829