



Material description: SULPHURIC ACID 5.0M (10N)

CAS No: 7664-93-9

Supplied as:	Product code	Description
	CHE3680	SULPHURIC ACID 5.0M (10N) 2.5L.
	CHE3682	SULPHURIC ACID 5.0M (10N) 5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	mol/L	4.995	5.005		

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2543/MS/19890101