



Material description: SULPHURIC ACID 90% w/w (for milk testing)  
CAS No: 7664-93-9

Supplied as: Product code Description  
CHE3708 SULPHURIC ACID 90% w/w (for milk testing) 2.5L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Strength	% w/w	89.5	91.5		
Wt. per ml.	g	1.812	1.818		
Chloride	%		0.00002		
Nitrate	%		0.00002		

Revision date: 25/03/2009

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 2400/MS/20090325