

Material description: 2,2,4-TRIMETHYLPENTANE pure
CAS No: 540-84-1

Supplied as: Product code Description
CHE3736 2,2,4-TRIMETHYLPENTANE pure 500ml.
CHE3738 2,2,4-TRIMETHYLPENTANE pure 2.5L.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Purity	%	99			
Boiling range	°C	97	99		
Wt. per ml.	g			0.69	
Absorbance @ 200nm	%		1.000		
Absorbance @ 212nm	%		0.300		
Absorbance @ 220nm	%		0.150		
Absorbance @ 240nm	%		0.040		
Absorbance @ 250nm	%		0.020		
Absorbance @ 260nm	%		0.015		

Revision date: 11/11/2020 Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1656/MS/20201111