

Material description: TALC  
CAS No: 14807-96-6

*Supplied as:* Product code            Description  
                         CHE3740                    TALC 1Kg.  
                         CHE3742                    TALC 3Kg.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Magnesium Oxide; MgO	%			29.7	
Silicon Dioxide; SiO <sub>2</sub>	%			52	
Aluminium Oxide; Al <sub>2</sub> O <sub>3</sub>	%			1.1	
Calcium Oxide; CaO	%			4.8	
Loss on ignition	%			11.4	
Acid soluble material	%			16	
Solubility in water	%			0.1	
pH	pH	7.0	9.0		
Asbestos type materials					Not detected

Revision date: 19/03/2009                    Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1815/MS/20090319