

Material description: TETRACHLOROETHYLENE pure  
CAS No: 127-18-4

<i>Supplied as:</i>	Product code	Description
	CHE3752	TETRACHLOROETHYLENE pure 500ml.
	CHE3754	TETRACHLOROETHYLENE pure 2.5L.
	CHE3756	TETRACHLOROETHYLENE pure 25L.

## Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Boiling range	°C	120	122		
Wt. per ml.	g	1.618	1.622		
Non-volatiles	%		0.001		

Revision date: 11/08/2017

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 0658/MS/20170811