



**SCIENTIFIC
LABORATORY
SUPPLIES**

Material description: THYMOL BLUE

CAS No: 76-61-9

Supplied as: Product code Description
CHE3778 THYMOL BLUE 5g.

Specification

Parameter	Units	Minimum value	Maximum value	Typical value	Notes
Colour change; Red to Yellow					pH 1.2-2.8
Colour change; Yellow to Blue					ph 8.0-9.6

Revision date: 01/01/1989

Issue date: 17/09/2024

Specification correct at time of issue but is subject to change without notice. It is the users responsibility to verify the material is suitable for the required application.

Ref: 1861/MS/19890101